

Children's Hospital Los Angeles  
 Alexander R. Judkins, MD  
 Department of Pathology & Laboratory Medicine  
 Pathologist-in-Chief and Laboratory Director  
 Phone: 323.361.2469, 323.361.2426  
 Fax: 323.361.8004  
 CLIA Number: 05D0542989  
 California State License CLF260  
 CAP Number: 2266301



**Ship To:**  
 Department of Pathology and Laboratory Medicine  
 Children's Hospital Los Angeles  
 6450 Sunset Blvd.  
 MRI/Weingart, 2nd Floor, Room 2-11  
 Los Angeles, CA 90027

**IMMUNOHISTOCHEMISTRY, SPECIAL STAINS & ELECTRON MICROSCOPY REQUISITION**

All information must be completed before sample can be processed.

**PATIENT INFORMATION**

\_\_\_\_\_  
 Last Name First Name MI

DOB (MM/DD/YYYY): \_\_\_\_\_ Gender:  M  F  Unknown

MRN: \_\_\_\_\_

**BILLING INFORMATION**

Referring Institution \_\_\_\_\_

CHLA Account Number:\* \_\_\_\_\_

Hospital/Laboratory Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip Code: \_\_\_\_\_

Accounts Payable Contact Name: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

**CLINICAL AND SPECIMEN INFORMATION**

**Services Requested**

- Stain Only
- Stain with Interpretation

**Number of Samples Submitted:**

- Paraffin Block(s): \_\_\_\_\_
- Slides (two per stain requested): \_\_\_\_\_
- Frozen Tissue(s) (OCT, Snap Frozen): \_\_\_\_\_

- Check here if you would like unstained slides returned

**Fixative:**

- Formalin
- B5
- Glutaraldehyde
- Other: \_\_\_\_\_

**Accession Number:** \_\_\_\_\_

**Block Number (s):\*** \_\_\_\_\_

**SEE PAGE 2 FOR SAMPLE REQUIREMENTS AND SHIPPING INSTRUCTIONS.**

\*See reverse side to open an account with CHLA Laboratory.

**IHC, SPECIAL STAINS & ELECTRON MICROSCOPY MENU**

- |   |  |   |
|---|--|---|
| <ul style="list-style-type: none"> <li><input type="checkbox"/> ABCBII (BSEP)</li> <li><input type="checkbox"/> ADV</li> <li><input type="checkbox"/> AFP</li> <li><input type="checkbox"/> ALK</li> <li><input type="checkbox"/> AP-2 beta</li> <li><input type="checkbox"/> BCL2</li> <li><input type="checkbox"/> BCL6</li> <li><input type="checkbox"/> BCOR</li> <li><input type="checkbox"/> Beta-Catenin</li> <li><input type="checkbox"/> Brachyury</li> <li><input type="checkbox"/> BRG1 (SMARCA4)</li> <li><input type="checkbox"/> C4d (paraffin)</li> <li><input type="checkbox"/> Calcitonin</li> <li><input type="checkbox"/> Calponin</li> <li><input type="checkbox"/> Calretinin</li> <li><input type="checkbox"/> CD1a</li> <li><input type="checkbox"/> CD2</li> <li><input type="checkbox"/> CD3</li> <li><input type="checkbox"/> CD4</li> <li><input type="checkbox"/> CD5</li> <li><input type="checkbox"/> CD7</li> <li><input type="checkbox"/> CD8</li> <li><input type="checkbox"/> CD10</li> <li><input type="checkbox"/> CD15</li> <li><input type="checkbox"/> CD19</li> <li><input type="checkbox"/> CD20</li> <li><input type="checkbox"/> CD21</li> <li><input type="checkbox"/> CD30</li> <li><input type="checkbox"/> CD31</li> <li><input type="checkbox"/> CD34</li> <li><input type="checkbox"/> CD43</li> </ul> | <ul style="list-style-type: none"> <li><input type="checkbox"/> CD45/LCA</li> <li><input type="checkbox"/> CD56</li> <li><input type="checkbox"/> CD61</li> <li><input type="checkbox"/> CD68/KPI</li> <li><input type="checkbox"/> CD79a</li> <li><input type="checkbox"/> CD99/MIC2</li> <li><input type="checkbox"/> CD117</li> <li><input type="checkbox"/> CD138/Syndecan</li> <li><input type="checkbox"/> CD163</li> <li><input type="checkbox"/> Chromogranin A</li> <li><input type="checkbox"/> CK (AE1/AE3)</li> <li><input type="checkbox"/> CK (CAM 5.2)</li> <li><input type="checkbox"/> CK7</li> <li><input type="checkbox"/> CK19</li> <li><input type="checkbox"/> CK20</li> <li><input type="checkbox"/> CMV</li> <li><input type="checkbox"/> C-Myc</li> <li><input type="checkbox"/> Collagen IV</li> <li><input type="checkbox"/> D2-40/Podoplanin</li> <li><input type="checkbox"/> Desmin</li> <li><input type="checkbox"/> E-Cadherin</li> <li><input type="checkbox"/> EBER (ISH)</li> <li><input type="checkbox"/> EBV:(LMP-1)</li> <li><input type="checkbox"/> EMA</li> <li><input type="checkbox"/> Factor XIIIa</li> <li><input type="checkbox"/> Fascin</li> <li><input type="checkbox"/> FLI-1</li> <li><input type="checkbox"/> GFAP</li> <li><input type="checkbox"/> GLUT-1</li> <li><input type="checkbox"/> Glutamine Synthetase</li> <li><input type="checkbox"/> Glypican-3</li> </ul> | <ul style="list-style-type: none"> <li><input type="checkbox"/> H3K27M</li> <li><input type="checkbox"/> H3K27me3</li> <li><input type="checkbox"/> HCG</li> <li><input type="checkbox"/> Histone H3.3 G34R</li> <li><input type="checkbox"/> HMB45</li> <li><input type="checkbox"/> HMGA2</li> <li><input type="checkbox"/> Helicobacter pylori</li> <li><input type="checkbox"/> Hep-Par1</li> <li><input type="checkbox"/> HSV1</li> <li><input type="checkbox"/> HSV2</li> <li><input type="checkbox"/> IDH1</li> <li><input type="checkbox"/> Inhibin Alpha</li> <li><input type="checkbox"/> INI-1</li> <li><input type="checkbox"/> Kappa (ISH)</li> <li><input type="checkbox"/> Ki67</li> <li><input type="checkbox"/> L1CAM</li> <li><input type="checkbox"/> Langerin</li> <li><input type="checkbox"/> LIN28</li> <li><input type="checkbox"/> LEF1</li> <li><input type="checkbox"/> Lysozyme/Muramidase</li> <li><input type="checkbox"/> MART -1/Melan A</li> <li><input type="checkbox"/> MDR3</li> <li><input type="checkbox"/> MGMT</li> <li><input type="checkbox"/> Myeloperoxidase</li> <li><input type="checkbox"/> MSA</li> <li><input type="checkbox"/> MUC4</li> <li><input type="checkbox"/> MUM1</li> <li><input type="checkbox"/> Myf-4/Myogenin</li> <li><input type="checkbox"/> MyoD1</li> <li><input type="checkbox"/> Neurofilament 200kD</li> <li><input type="checkbox"/> Neurofilament (2F11)</li> </ul> |
|---|--|---|

**Immunofluorescent Stains (frozen sections required)**

- C1q
- C3c
- C4d
- Fibrinogen
- IgA
- IgG
- IgM

**IHC Stains for Muscle Biopsies**

- Alpha Sarcoglycan
- Beta-Dystroglycan
- Dystroglycan (C-Terminus)
- Dystroglycan (Rod- Domain)
- Dystroglycan (N-Terminus)
- Dysferlin
- Merosin
- Spectrin
- Myosin Heavy Chain (Fast)
- Myosin Heavy Chain (Slow)
- Myosin Heavy Chain (Developmental)
- Myosin Heavy Chain (Neonatal)

**Muscle Biopsy Enzyme Histochemistry/Special Stains**

- Gomori's One-Step Trichrome
- Oil Red O
- PAS
- PAS Diastase
- Acid Phosphatase
- Alkaline Phosphatase
- ATPase (4.3, 4.6, 9.4)
- Cytochrome Oxidase
- NADH
- SDH

**Special Stains**

- AFB
- Alcian Blue 2.5
- Congo Red
- Copper
- Elastic Van Gieson
- Fite
- Fontana-Masson
- Giemsa-May Grünwald
- GMS
- Gram
- Iron
- Jones
- Luxol Fast Blue
- Mucicarmine
- Oil Red O
- PAS
- PAS Diastase
- Reticulin
- Warthin-Starry
- Trichrome (Masson's)
- Von Kossa

**Electron Microscopy**

- Full EM with interpretation
- Full EM without interpretation (Images provided to requesting Institution without interpretation)
- Process Only (create Epon block)
- Process+Semithin Sections
- Scoping only (grids provided by submitting institution)



**Ship To:**

Department of Pathology and Laboratory Medicine  
Children's Hospital Los Angeles  
4650 Sunset Blvd. MS#43  
MRI/Weingart, 2nd Floor, Room 2-11  
Los Angeles, CA 90027

**TEST REQUISITION REQUIREMENTS**

1. All samples should be clearly labeled with **patients first and last name, date of birth, hospital/ laboratory name, and block or slide number.** Failure to fully complete the test requisition may delay analysis.
2. **All fields with red text are required.**
3. We will notify you within 24 hours of receipt if we are unable to perform testing due to failed specimen integrity. Any specimen determined to be inadequate for testing will be discarded after 48 hours.
4. **Cancellations:**

**SHIPPING AND HANDLING INSTRUCTIONS**

**FRESH FROZEN TISSUE:**

1. Ship fresh frozen tissue on dry ice with a minimum of 5 kg of dry ice. Overnight delivery required to prevent thawing of tissue.

**PARAFFIN BLOCK(S) OR UNSTAINED SLIDES :**

1. Ship paraffin blocks or slides at room temperature.

**GENERAL INSTRUCTIONS:**

1. We will notify you within 24 hours of receipt if we are unable to perform testing due to compromised sample integrity.
2. Please notify us ASAP in writing if you wish to cancel a test. Cancellations cannot be accepted once testing has been initiated.
3. We accept samples Monday through Thursday from 7:00 AM to 4:00 PM PST. We also accept samples on Friday by 11:00 AM PST. All packages should be mailed for receipt by Friday. Holidays and weekends should be taken into consideration before mailing samples.
4. To ensure sample integrity, use of the following delivery priorities is highly recommended.  
FedEx: First Overnight  
UPS: Next Day Air Early AM
5. **Your specimen is important to us. Please email the tracking number to [PLMTrack@chla.usc.edu](mailto:PLMTrack@chla.usc.edu) at the time of shipment and include contact information to be used in the event your sample is not received.**

**BILLING INFORMATION**

1. For billing inquiries, please call (877) 543-9522.
2. If you are interested in opening an account with Children's Hospital Los Angeles, please contact our Laboratory Service Center at (877)543-9522. Please be prepared to provide the following information:
  - a. Name of Institution
  - b. Address
  - c. Phone/Fax Number
  - d. Laboratory Contact Name and phone number
  - e. Accounts Payable Contact Name and phone number
3. Third party billing is not offered at this time.

**CHILDREN'S CONNECT**

Children's Connect is a web-based portal providing 24/7 access to laboratory test order entry and results retrieval.

To request access or to receive more information, please contact us at:  
(877) 543-9522



**CONTACT US**

For all other inquiries, please contact our Laboratory Service Center at:  
**(877)KIDZ-LAB or (877) 543-9522**  
or visit our website at:  
**CHLA.org/CPM**